



InVivo Biosystems

Right Model. Right Test. Right Insights.

InVivo Longevity Platform

EXAMPLE SERVICE REPORT



Customer	ABC Aging
Project contact	John Doe
Project code	ABCA-001
Project scope	<ol style="list-style-type: none">1. Understand whether compound Lu0128 has an effect on lifespan and healthspan.2. Identify all the genes whose expression is affected by compound Lu0128.3. Identify pathways that are targeted by compound Lu01284. Determine the extent to which known aging-related pathways are targeted by compound Lu0128.



Table of contents

Report summary	3
Project Scope and Aims	3
Executive Summary	3
Project Conclusions	4
Results / Analyzed Data	5
Step I. Dosage, toxicity, and feasibility.	5
1.1 Stress response assay and dose selection	5
1.2 Viability, Growth, and Development assay to determine toxicity in <i>C. elegans</i> .	6
Step II. Does compound Lu0128 extend lifespan and/or healthspan?	7
2.1. Lifespan Assay	7
2.2. Healthspan, movement, and morphology analysis	11
Step III. Whole Transcriptome analysis	14
3.1 Gene expression analysis	14
3.2 Mechanism of Action report	16
Methods, Supplemental Materials, and Raw Data	17

Report summary

Project Scope and Aims

Use the nematode *C. elegans* as a model to gain insight into the impact of compounds Lu0128 on lifespan, healthspan, gene expression and aging-specific pathways.

Using a 3-pronged approach, we:

1. Measured toxicity and determined dosage of Lu0128 for lifespan/healthspan studies.
2. Measured the effect of compounds Lu0128 on lifespan.
3. Identified Lu0128's mode of action.

Step 1: We determined a physiologically effective dosage and delivery method for Lu0128, and assessed the toxicity of Lu0128 to *C. elegans*.

Step 2: We determined whether Lu0128 extends lifespan and healthspan using automated survival measurement and in-depth healthspan analysis.

Step 3: We Identified genes and pathways related to Lu0128-mediated pro-longevity effects using whole transcription analysis.

Executive Summary

Oxidative stress has been shown to play a role in aging and age-related illness. ABC Aging's product Lu0128 is a derivative of Glutathione that has been modified for enhanced potency and biostability. The nematode worm *C. elegans* is an established and validated model system for testing products that modulate lifespan or treat age-related illness. Lu0128 significantly increased the average lifespan of *C. elegans*, similar to the well-established life-extending compound that we used as a positive control. Furthermore, deep phenotyping revealed that Lu0128 improved healthspan relative to both vehicle (no treatment) and positive controls. Whole transcription analysis indicates that the effect of Lu0128 on lifespan and healthspan likely results from modulation of oxidative stress, consistent with this class of compounds.

Project Conclusions

Step I. Dosage, Toxicity, and Feasibility tests

- Compound Lu0128 is water-soluble and can be delivered to worms by permeating solid agar media with an aqueous solution.
- Based on the dose response curves of two stress biosensors, a working EC_{50} of $21\mu\text{M}$ (rounded to $20\mu\text{M}$) was selected as a benchmark for further testing.
- Growth, viability, and toxicity assays indicated a maximum non-toxic dosage of $10\mu\text{M}$ Lu0128.
- Recommend proceeding with a dosage of $10\mu\text{M}$ Lu0128 for lifespan assay.

Step II. Lifespan assay and survival analysis.

- Worms treated with Lu0128 showed an increased median and maximum lifespan
- Worms treated with $10\mu\text{M}$ Lu0128 showed significantly increased lifespan relative to vehicle control (Log-rank test, P-value = <0.0001).
- Worms treated with Lu0128 remained in a mobile state longer than those treated with vehicle control.
- Principal components analysis revealed that the movement features of worms treated with Lu0128 clustered exclusively from those treated with vehicle control and similarly to those treated with the positive control.

Step III. Whole transcriptome analysis

- A total of 74 genes showed significantly ($p < 0.00001$) different expression at day 10 of adulthood in Lu0128 treated vs. untreated animals.
- Lu0128 showed a significant change in the activity of pathways linked to oxidative phosphorylation and oxidative stress.

Summary conclusions for this project

1. Lu0128 significantly increased lifespan and healthspan.
2. Lu0128 possibly influences lifespan through modulation of oxidative phosphorylation and oxidative stress.



Results / Analyzed Data

Step I. Dosage, toxicity, and feasibility.


The goal of this set of experiments is to determine experimental conditions for steps II and III.

1.1 Stress response assay and dose selection

The EC_{50} is the concentration of a drug that gives half-maximal response. We determine a working EC_{50} for Lu0128 by measuring physiological stress in response to varying concentrations using our fluorescent stress biosensor *C. elegans* strains.

Transgenic biosensor worms express red fluorescent protein (RFP) in response to stress and constitutively express green fluorescent protein (GFP).

- Day 1 adult biosensor worms were treated with a 1 μ M-10 mM range of Lu0128 concentrations.
- Fluorescent images were captured after 24 hours of treatment.
- Stress response was quantified by normalizing the RFP signal to the GFP signal.



Compound Lu0128 is water-soluble and can be delivered to worms by permeating solid agar media with an aqueous solution. Based on the dose response curves of two stress biosensors (Figure 1.1.), a working EC_{50} of 21 μ M (rounded to 20 μ M) was selected for further testing.

1.2 Viability, Growth, and Development assay to determine toxicity in *C. elegans*.

- Worms were exposed to three concentrations of compounds (EC_{25} , EC_{50} , EC_{100}) from embryo to first day of adulthood.
- 50 treated embryos were plated on each of three replicate dishes.
- Quantity and size of worms was measured each day through day 1 of adulthood.




Figure 1.2. Viability, Growth, and Development assay to determine toxicity of Lu0128 in *C. elegans*. A) Violin plot showing distribution of worm lengths during development from larvae to adult when treated with varying concentrations of Lu0128. B) Day three snapshot of worm development including embryonic lethality (black). EC₂₅, EC₅₀, EC₁₀₀ are based on dose-response curves in Experiment 1.1 and correspond to concentrations of 10µM, 20µM, and 40µM.

In viability assays, Lu0128 exposure at the EC₁₀₀ caused a clear developmental delay (Figure 1.2.A.) and reduced viability (Figure 1.2.B.) . A slight developmental delay was observed at the EC₅₀. Worms treated with 10µM Lu0128 (EC₂₅ = 0.5 x EC₅₀) showed normal viability and development. We therefore recommended using a dosage of 10µM Lu0128 for Step II. This recommendation was approved by ABC Aging.


Step II. Does compound Lu0128 extend lifespan and/or healthspan?

2.1. Lifespan Assay

- Conditions tested:
 - Compound Lu0128 at 10µM
 - N-acetyl-L-cysteine (NAC), a class-specific positive control, at 50µM
 - Vehicle control: water

- 3 biological replicates of 150 worms each derived from independent synchronizations. >450 worms total.
- 10 μ M Lu0128 delivered by permeating solid agar media and food.
- Treatment with Lu0128 began on day 1 of adulthood (day 3 post-hatch)
- Automated lifespan recording was performed using Life and Death Instrument (LaDI, InVivo Biosystems, Inc)

To understand the physiological action and therapeutic potential of compounds targeting aging, it is crucial to use quantitative models. The analysis of lifespan data is grounded in the study of two mathematical functions: the survival curve and the hazard function. The survival curve (Figure 2.1.A.) describes the fraction of a tested population that remains alive over time. The hazard function (Figure 2.1.B.) is related to the survival curve and provides an intuitive measure of the risk of death. This function describes the probability that a typical individual who is currently alive will soon die, providing a clear visualization of the way a compound may change patterns in mortality.



Worms that receive the Lu0128 treatment (blue line) lived significantly longer (maximum lifespan 33 days) than worms that did not (green line, maximum lifespan 28 days). Treatment with positive control NAC (orange line) also increased lifespan as expected. Tables 2.1 to 2.4 detail the raw measurements, calculations and statistical analysis related to the survival curve in Figure 2.1.A.

Lu0128 treatment temporally shifted survival by delaying early deaths (Figure 2.1.A). This is reflected in a delayed increase in the hazard rate (Figure 2.2.B.). Lu0128 increased maximum lifespan (Table 2.1.), with a small subpopulation exceeding the maximum lifespan of the

positive control. The hazard ratio (Table 2.4.) compares two treatments. It is not computed at any one time point, but includes all the data in the survival curve. The hazard ratio between the Lu0128-treated population and the untreated control population is 0.46, indicating that the rate of deaths in the Lu0128-treated group is 54% lower than in the control group.

	Lu0128 (10µM)	Vehicle (H ₂ O)	Positive control
Death events	411	431	428
Censored subjects (#worms)	N/A	N/A	N/A
Median lifespan (days)	24	21	25
Average lifespan (days)	22	19	23
Maximum lifespan (days)	33	28	35

treatment	Mean (days)	Mean 95 % CI	25% mortality	50%	75%	90%	100%
Lu0128	21.49	20.93 - 22.06	19.03	20.04	21.25	24.14	33
Vehicle (water)	19.81	19.43 - 20.19	17.18	19.16	20.31	20.99	28
NAC	22.46	21.98 - 22.95	18.5	21.57	23.19	25.08	35

Mean is calculated from the area under the survival curve. Age at percent mortality is calculated from a linear interpolation of the survival curve.

Treatment	at 25% mortality	at 50%	at 75%	at 90%
Lu0128 vs. Vehicle	0.002 **	3.90E-7 ****	0.0042 **	0.0052 *
Lu0128 vs. NAC	0.6782	0.0293 *	0.0593	0.4316 *
NAC v.s. Vehicle	0.0001***	0.00236 **	0.0044 **	0.00454 **

p-values for Fisher's exact test calculating the probability of obtaining the observed data at different stages of mortality if there is no difference between the two populations. P-values shown.

Curve comparison	Hazard Ratio	Log-rank (Mantel-Cox) test [†]	Gehan-Breslow-Wilcoxon test [‡]
Lu0128 (10µM) vs vehicle.	0.46	<0.0001 ****	<0.0001 ****
Lu0128 (10µM) vs positive control	1.07	0.17	0.23

The hazard ratio represents the overall rate of death relative to control. It is generalized to all of the data in the survival curve, and will not be consistent over time for survival curves that cross or inflect.
[†]Mantel-Cox test compares groups across the duration of the lifespan.
[‡]Gehan-Breslow-Wilcoxon test applies more weight to earlier deaths.
 Numbers and asterisks represent P-value and significance, respectively.

Worms treated with 10µM Lu0128 showed

- An increased median and maximum lifespan
- A significantly increased lifespan relative to vehicle control (Log-rank test, P-value = <0.0001).
- A significantly lower hazard ratio than untreated worms.

2.2. Healthspan, movement, and morphology analysis

To measure healthspan, precise movement data was collected concurrently with survival data from lifespan assay. Images collected automatically from the LaDI for lifespan recording were analyzed to assess healthspan.

- Aggregate movement data was captured from LaDI analysis

Worms treated with Lu0128 remained in a mobile state significantly longer than those treated with vehicle control and slightly longer than those treated with positive control.




Figure 2.2. Aggregate movement analysis over duration of lifespan. Black line: Survival
Green line: Aggregate mean displacement of worms on the plate. A single trace from one condition is shown for demonstration.

Treatment	HS ₅₀ P-value	Log-rank test
Lu0128 (10µM) vs vehicle	<0.0001 ****	<0.0001 ****
Lu0128 (10µM) vs positive control	0.013 *	0.021 *
Positive control vs vehicle	<0.0001 ****	<0.0001 ****

Healthspan-50 (HS₅₀) is the time (days) at which aggregate movement is reduced to 50% of maximum. HS₅₀ = 12.5, 14.5, and 15.0 for vehicle, positive control, and Lu0128, respectively.
Lox-rank test compares the aggregate movement curves over the duration of the lifespan.

- Deep phenotyping: Worm movement and morphology features were extracted from the images and subjected to principal components analysis.

The automation of the lifespan assay allows the simultaneously gathering of over 20 different variables related to the health of the observed populations. In order to reduce the dimensionality of such a rich data set, we perform a Principal Component Analysis (PCA), a technique that allows the compression and classification of data to find a new, smaller set of variables, that nonetheless retains most of the sample's information. We apply PCA to extract and compare 7 measurements of Principal components (PCs) of vitality between the populations. Such measurements include morphology and movement. Here, the PCA revealed that the movement features of worms treated with Lu0128 clustered exclusively from those treated with vehicle control and similarly to those treated with the positive control NAC (Figure 2.4.).




Figure 2.4. Principal components analysis of healthspan and movement features. Each dot represents one worm. A random subset of twenty worms from each treatment is shown for clarity. Features include worm area, length, width, curvature, aggregate movement, among others. Yellow: vehicle control. Red: positive control. Blue: Lu0128.

Step III. Whole Transcriptome analysis

To determine the mode of action for Lu0128-mediated lifespan extension, we analyzed gene transcription in young and aged adults by RNA-Seq (whole transcription analysis). For each condition, >300 worms were treated in parallel with the lifespan assay (Step II) and maintained under identical conditions.

- Experimental conditions:
 1. Lu0128 10 μ M Day 1 adult
 2. Lu0128 10 μ M Day 15 adult
 3. Vehicle Day 1 adult
 4. Vehicle Day 15 adult

3.1 Gene expression analysis

A total of 74 genes showed significantly ($p < 0.00001$) different expressions at day 10 of adulthood in Lu0128 treated vs. untreated animals (Figure 3.1, Table 3.1. and Supplementary Materials).





Figure 3.1. Volcano plot shows the negative log₁₀ of P-values as a function of log₂ fold changes of *C. elegans* genes for Lu0128 treated vs. untreated day 10 adult animals. Lu0128 treated 10 day old animals have gene expression changes compared with untreated 10 day old animals. Red dots indicate genes with significant changes.

Table 3.1 Individual genes highly up- or downregulated on day 10 by Lu0128. Top 20 most significant hits shown: Comprehensive list of genes is enclosed in Supplemental Materials.

gene	logFC	logCPM	LR	PValue
Y73B3A.1	-9.25	5.94	256.75	8.76E-58
asp-12	-6.87	8.85	164.32	1.29E-37
cTel55X.1	7.25	3.82	94.51	2.44E-22
cyp-35A2	-2.78	5.65	91.62	1.05E-21
mtl-2	-3.16	5.77	77.76	1.17E-18
Y73B3A.3	-10.58	2.53	71.12	3.36E-17
lea-1	2.69	9.58	67.21	2.43E-16
Y73B3A.21	-10.31	2.23	65.09	7.16E-16
unc-18	-3.10	4.08	61.51	4.41E-15
ugt-22	-1.81	6.47	60.93	5.93E-15
CE7X_3.1	-9.78	1.80	51.18	8.44E-13
T28F3.8	-1.99	5.41	50.78	1.03E-12
Y73B3A.20	-9.50	1.29	45.46	1.56E-11
F31F7.1	-1.75	6.37	45.32	1.67E-11
F23A7.8	-1.73	5.74	43.06	5.30E-11
F42A9.6	1.62	5.66	41.28	1.32E-10
cyp-35A3	-4.98	2.24	40.23	2.25E-10
cth-1	-2.36	6.33	39.57	3.17E-10
atic-1	-1.81	5.42	38.14	6.58E-10

3.2 Mechanism of Action report

Our pathway analysis shows that worms treated with Lu0128 exhibit a significant change in the activity of pathways linked to oxidative phosphorylation and oxidative stress.



Methods, Supplemental Materials, and Raw Data

The following files were sent to the customer via email in a .zip folder

1. Assay design, materials and methods
 - ABCA-001_Materials_Methods.pdf
2. Raw data for feasibility and dosage tests
 - ABCA-001_doseage_avoidance.xlsx
3. Raw lifespan/healthspan assay data in spreadsheet.
 - ABCA-001_lifespan_scoring.csv
 - ABCA-001_movement.csv
4. Raw RNA-seq data
 - ABCA-001_day1.csv
 - ABCA-001_day8.scv
5. Supplemental data
 - ABCA-001_supplement.pdf