



## InVivo FLite Dinolite Microscope Working Distance and Field of View

<b>Working distance and field of view (inches)</b>														
Magnification	10	20	30	40	50	60	70	80	90	100	110	120	130	140
Working Distance	9.354	4.528	2.969	2.232	1.827	1.579	1.429	1.339	1.283	1.256	1.248	1.256	1.276	1.307
Field of View	1.539	0.768	0.512	0.386	0.307	0.256	0.220	0.193	0.169	0.154	0.142	0.130	0.118	0.110
	x 1.142	x 0.579	x 0.386	x 0.287	x 0.232	x 0.193	x 0.165	x 0.146	x 0.130	x 0.114	x 0.106	x 0.094	x 0.091	x 0.083
Depth of Field	0.969	0.263	0.126	0.076	0.052	0.039	0.03	0.024	0.02	0.017	0.015	0.013	0.012	0.011