

InVivo PhysioTemp Benchtop Incubator



- 1 Power Button** turns incubator on or off
- 2 Temperature Indicator** reports "set" temperature & "chamber" temperature
- 3 Temperature Adjustment** is in increments of 1C & changed by pressing "up" or "down"
- 4 Chamber Light Button** turns the interior LED ON or OFF
- 5 Door Latch** opens and closes the door
- 6 Passthrough** used for power or USB cables

Specifications

Temperature Range: 5C-50C

Temperature Accuracy: 0.5C

Temperature Uniformity (shelves): 0.5C

Temperature Uniformity (Chamber at drip tray): 1C@20C / 1.5C@37C



Temperature Increment: 1C

Electrical Input: 1C@20C / 1.5C@37C

Electrical Adapter: 120/220VAC - 1.5A w 12VDC Port Option

QC report

QC Function Tests					
Power On (Test 1)	PASS	LED On (Test 2)	PASS	QC Runtime	192 Hours
Heats (Test 3)	PASS	Cool (Test 4)	PASS	Min/Max	4.4c/51.2c

QC Stability Tests	
15c	<p>2 Days @ 15c</p> 
20c	<p>5 Days @ 20c</p> 
37c	<p>1 Day @ 37c</p> 