

# Rapid genotyping of live zebrafish embryos using the Zebrafish Embryo Genotyper device (ZEG)

Andy Willaert, PhD  
Center of Medical Genetics  
University Hospital Ghent  
Belgium



CENTRUM MEDISCHE  
GENETICA GENT



Zebrafish Facility Management  
Workshop, Lisbon, 2018



# Ghent Zebrafish Facility



- 6 “semi-closed” recirculating systems (ZebTEC and WTU systems, Tecniplast)
- Housing of 17,000 zebrafish
- ‘Danio Data’ (Fulcrum) zebrafish husbandry management software
- Rotifer-polyculture (Best *et al.*, 2010), Micro-artemia (Ocean Nutrition) and dry food (Skretting)



# Adult Zebrafish Genotyping

## Finclipping of adult zebrafish


Raise fish to adult age



Manual fin clip



DNA extraction/Genotyping




- Requires **space, effort and expenses** of raising more adult fish than needed to ensure obtaining sufficient numbers of animals of the desired genotype
- Finclipping procedure is **time- and labor intensive**
- **Early identification of genotypes** needed in experiments dealing with larval phenotypes and early death or in experiments that require sample pooling




Kosuta *et al.*, 2018  
Wilkinson *et al.*, 2013

## Larval fin clipping


Intact larval fin at 3dpf



Tail fin clipping at 3dp





Fin regrowth at 5dpf



- Sufficient amounts of DNA for **accurate genotyping**, **good survival rates**
- Several **critical steps within the protocol** must be correctly followed to obtain qualitative results (high-precision fin clip; tissue transfer to test tube must be assured)
- Finclipping procedure is **time- and labor intensive on a larger scale**



## Zebrafish Embryo Genotyper (ZEG)




*Lambert et al., 2018 (PlosOne)*







# Larval Zebrafish Genotyping

## Genetic material extraction using the ZEG



## Genetic testing of DNA obtained from the ZEG



- No apparent effects on body morphology and motor behavior at 7dpf or long-term growth and survival (>90%)



# ZEG testing

Success rate of PCR/genotyping using ZEG solution is high; some general conclusions/trends can be drawn:

- During the ZEG procedure, the embryos are unharmed, with no apparent adverse effects on survival (survival is usually 100%)
- The mean PCR/genotyping success rate is 96%. Only when no PCR amplicon can be detected on gel/fragment analyzer then sequencing fails





- ZEG works equally well for 2dpf and 3dpf embryos, with similar PCR/genotyping success rates. We now mostly use 2dpf embryos and obtain genotypes when fish are 4dpf



# ZEG testing

- 5µl of ZEG solution in a 20µL volume PCR works best (highest PCR/genotyping success rates) although 4 or 3 µL also works fine in most cases. Regular touch down PCR programs work well
- With increasing vibration times, more DNA is extracted from embryos



- High correlation between ZEG genotypes and whole embryo genotypes
- ZEG solution works equally well for genotyping with traditional Sanger sequencing and NGS-Miseq sequencing (eg to characterize mosaic embryos from CRISPR experiments)



- The ZEG device is an automated system that rapidly obtains DNA from zebrafish embryos for genotyping, while keeping the animals alive and intact
- Allows researchers to identify and raise mutants of interest at an early stage of development and therefore saves time, effort and money
- High PCR/sequencing success rate from 2dpf and 3dpf embryos with no or little harm to the embryos up to 30 minutes of vibration
- High correlation between ZEG genotypes and whole embryo genotypes



# Acknowledgements

**Center for Medical Genetics Ghent, Ghent  
University, Ghent, Belgium**

- Hanna De Saffel
- Annekatrien Boel
- Patrick Sips
- Paul Coucke



CENTRUM MEDISCHE  
GENETICA GENT



**Department of Pediatrics and Department of  
Engineering, University of Utah, Salt Lake City,  
Utah, United States of America**

- Joshua L. Bonkowsky
- Christopher J. Lambert



Startup Company 'Fluidx' (<https://www.wfluidx.com/>)



RESEARCH ARTICLE

An automated system for rapid cellular  
extraction from live zebrafish embryos and  
larvae: Development and application to  
genotyping

Christopher J. Lambert<sup>1</sup>, Briana C. Freshner<sup>2</sup>, Arlen Chung<sup>1</sup>, Tamara J. Stevenson<sup>2</sup>, D.  
Miranda Bowles<sup>2</sup>, Raheel Samuel<sup>1,3†</sup>, Bruce K. Gale<sup>1,3†</sup>, Joshua L. Bonkowsky<sup>2‡\*</sup>