




Wmicrotracker Specifications

- ✓ Dimensions: 22cm x 28cm x 9.1cm (8.66in x 11in x 3.6in)
- ✓ Power source: 9V DC, 1.5 Amp switching Power
- ✓ Connection: USB-B cable
- ✓ Microplate format adapters: 96-w. fl at bottom, 96 w. "U" shaped bottom
- ✓ Optimal operating environment: 10°C to 37°C, low humidity (<50%)

PACKING LIST:

The system includes the following accessories:



USB cable

Microplate format adapters:

- 96 w. flat bottom
- 96 w. "U" shaped bottom

Acquisition Software
available from:

www.nemamatrix.com

Microplate reader system

SYSTEM INSTALLATION MINIMUM REQUIREMENTS:

- ✓ PC compatible with the following minimum requirements:
- ✓ Pentium II processor or above (>1GHz clock)
- ✓ 512Mb of RAM memory
- ✓ 1 USB port
- ✓ Windows XP 32 bits (or higher) operating system
- ✓ 200Mb of free HD space (>10Gb free HD space recommended for real time data saving)