









NemaMetrix NemaMetrix

Larval Stage 2 360 - 380 μm

Larval Stage 3/4 500 - 900 μm

Adult Day 1 1.0 mm

Adult Day 6 1.7 mm

arynx • gonad is appearing at the alf body posterior half body

• bean shaped vulva

• embryos and eggs present

• few embryos or eggs present

dark body color

ScreenChip size	SC10	SC20	SC30	SC40	SC60
Room temperature	(no food exposure, L1 arrest)	(Time post food exposure)		(days after first eggs formed in the gonade)	
18°C	12 hr	26 hr	48 - 50 hr	0 - 5 days	8+ days
20°C	12 hr	24 hr	46 - 48 hr	0 - 4.5 days	8+ days
22°C	12 hr	22 hr	44 - 46 hr	0 - 4 days	7+ days